



OptiMax Flexi PCR plate

96-well, white, non-skirted, cuttable, 0.1 ml

Low profile: height 15.6 – 16.6 mm

Cat. No.	Pack Size
31196	box of 25

Description:

OptiMax Flexi PCR plate with marker label on top is made of ultra-flexible material with sturdy tubes and can be easily cut in plate parts, multiples of 8- or 12-well strips and multiples of single tubes. Ideal to be combined with OptiClear flat 8-cap strip (cat.no. 31108). White wells are opaque, significantly reducing signal loss through tube walls and maximizing reflection. OptiMax Flexi PCR plate is designed for PCR and qPCR applications and fits most cycler models.

Advantages:

- ✓ Cuttable: can be easily cut with scissors
- ✓ Uniform, thin walls deliver maximum and consistent heat transfer
- ✓ Made of optimized blend of PP, designed for PCR and qPCR applications
- ✓ High sensitivity for low copy number targets
- ✓ Low evaporation: loss of PCR volume less than 3 %
- ✓ DNA, RNA, DNase, RNase, pyrogen and metal free
- ✓ Extra low adhesion: specially treated to reduce protein binding
- ✓ Autoclavable

Data supplied with products are correct and reliable to the best of our knowledge and belief. Under no circumstances will Selecta Biotech SE be liable for the consequential damage arising from the use of its products. Selecta Biotech SE is liable exclusively and restrictively to replace defective products. Selecta Biotech SE warrants its products to be free of defects in material and workmanship. Selecta Biotech SE will replace any products that are found to be defective, at no cost.

Warranty and Disclaimer: This product shall comply with its relevant specification and be fit for its stated purpose, but Selecta Biotech gives no other warranty and makes no representation as to description or quality. Any such warranty or representation is excluded, to the fullest extent permitted by law. In particular, but without limiting the foregoing, Selecta Biotech shall not be liable for the failure of the product to comply with its relevant specification where such failure arises as a result of: (i) customer negligence or because the customer failed to follow any of the applicable technical data or safety sheets, standard user materials, use guidelines or any other information provided by Selecta Biotech as to the storage, transportation, handling, use or maintenance of the products or other good practice regarding the same, or (ii) the customer altering the products in any way without the prior written consent from Selecta Biotech, or (iii) the products differing from the relevant specification as a result of changes made to ensure their compliance with applicable statutory or regulatory requirements.

Nothing shall limit or exclude Selecta Biotech's liability for death or personal injury caused by its negligence, fraud or fraudulent misrepresentation or any matter in respect of which it would be unlawful for Selecta Biotech to exclude or restrict liability. Without limiting the foregoing, Selecta Biotech shall under no circumstances whatever be liable to the customer, whether in contract, tort (including negligence), breach of statutory duty, or otherwise, for any loss of profit, or any indirect or consequential loss arising under or in connection with the products and Selecta Biotech's total liability to the customer in respect of all other losses arising under or in connection with the product, whether in contract, tort (including negligence), breach of statutory duty, or otherwise, shall in no circumstances exceed the price of the products supplied in respect of which the liability has arisen.